


# Laboration 11

DQF-COSY  
(**phase sensitive**)

KR


## Theory

## DQF-COSY



In order to resolve small J-couplings fine digital resolution is required, which significantly increases the experimental time. In general, the DQF-COSY experiment is recommended if a higher resolution is desired.

The DQF-COSY pulse sequence



The pulse  $p1$  must be set to the appropriate  $90^\circ$  time found. Note that the DQF-COSY experiment is sensitive to high pulse-repetition rate, i.e., it is important to choose a long recycle delay time  $d1$  in order to avoid multiple-quantum artifacts in the spectrum. A suitable value for this sample is  $d1 = 3$  sec.

## Practical

1. Run a *PROTON* according to "Bruker run manual for 500 MHz NMR".  
The 2-D exp should be recorded non spinning, so for a better result shim also the non spinning shims x, xz, y, yz
2. Optimize the values of **o1p** and **sw**
3. Check the  $^1\text{H}$  90°-pulse (lab 1).

## Experiment setup

1. **edc** or **new** and read the parameter set *COSYDQFPHSW* under Experiment. You will find it when you are in the correct directory. Press the arrow and choose the pathway to Brukers parametersets under Experiment Dirs. /opt/topspin21/exp/stan/nmr/**par**
2. **getprosol** (get probe and solvent specific parameters from prosol)

or if the  $^1\text{H}$  90°-pulse value need to be changed

Set the measured **p1** ( $^1\text{H}$  90°-pulse)

**getprosol 1H 'p1-value' 'p11-value'** (get probe and solvent specific parameters and use your adjusted p1 value to calculate related pulses)

3. If required, any acquisition parameter can be modified manually or in the *AcquPars* section, you can see what is valid for the parameters in *PulseProg*.



- a. **o1p** is the center of the 1H spectrum
- b. **sw** (F1) same exact value as sw (F2), you get it from 1H spectrum
- c. **td** (F2) is the time domain in the F2 dimension (usually set to 1K-2K)  
Set it to **2024** (2K)
- d. **td** (F1) is the number of experiments/increments to be recorded in the F1 dimension (usually set to 64w-512w)  
Set it to **512**
- e. Set appropriate **ns**, depending on the time you want the experiment to take.

4. turn of f the spinner
5. **rga**
6. **zg**

### **Process recorded data**

7. **xfb** add a window function and Fourier transform the data.
8. **abs1, abs2** perform a baseline correction
9. Phase, if necessary.
10. Set the reference. For external referencing go to the *Procpar*s section set sr(F1) and sr(F2) to the same value as in the <sup>1</sup>H spectrum
11. The resulting 2D spectrum can be also optionally symmetrized by using the **syma** command. Be careful! Why?